



Recombinant *Oryza sativa* subsp. japonica 12-oxophytodienoate reductase 7 (OPR7)

Product Code	CSB-BP765049OFG
Abbreviation	OPR7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6Z965
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MDRPPPDQQR QKQAPLFSPY QMPRFRLNHR VVLAPMTRCR AIGGVPGPAL AEYYAQRRTTQ GLLISEGTV VSPAGPGFPV VPGIYNQEQT DAWKKVVDVA HAKGGIFFCQ LWHVGRASHQ VYQPNGAAPI SSTDKPISAR WRILMPDGSY GKYPKPRRLA ASEIPEIVEQ YRQAAINAIE AGFDGIEIHG AHGYIIDQFL KDGINDRTDE YGGSLSNRCR FLLEVTRAVV SAIGADRVAV RISPAIDHLD AYDSDPIKLG MAVVERLNAL QQQSGRLAYL HVTQPRYTAY GQTESGQHGS AEEESRLMRT LRGTYQGTFM CSGGYTRELG LEAVESGDAD LVSYGRLFIS NPDLVERFRL NAGLNKYVRK TFYTPDPVVG YTDYPFLGQP KSRM
Source	Baculovirus
Target Names	OPR7
Protein Names	Recommended name: 12-oxophytodienoate reductase 7 EC= 1.3.1.42 Alternative name(s): 12-oxophytodienoate-10,11-reductase 7 Short name= OPDA-reductase 7 Short name= OsOPR7
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.