



Recombinant *Oryza sativa* subsp. japonica Probable gamma-aminobutyrate transaminase 4 (GABA-T)

Product Code	CSB-BP765052OFG
Abbreviation	GABA-T
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZH29
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MTTPHGLLQY SSGAFSDQVP ADDSAEEHGV KDHAMLAPFT AAWQTAISPP LVIERSEGCY VYDVNGTKYL DALAGLLSTA LGGSEPRLVK AATEQLNKLP FYHSFWNHTT RPSLDLAKEL ISMFTAREMG KVFFTNNGSE ANDSQVKIWW YYNNALGRP KKNIIISRTQS YHGTTFISAS LSGPLTLHQD FDLPGRFVLH TDCPHYWRFH LPGETEEFEFA TRLADNLENL ILKQGPETIA AFIAEPVIGA GGVILPPKTY FEKIQAVVKK YDILFIVDEV ITGFGR LGTM FGSDLYNIKP DLVSLAKALS SAYAPIGAIL VSPEISDVIH SHSNKLGTF A HGFTYSGHPV SCAVALEALK IYRERDIPGH VTHVAQR FQE GIKAF AAGSP IVGETR GVGL LIATEFTDNK SPYELFPFEW GVGEIFGQEC KKRGM MVKVL GNLIAMSPPL IITREEIDKL VSIYGEALKA TEERVAELKS KKN
Source	Baculovirus
Target Names	GABA-T
Protein Names	Recommended name: Probable gamma-aminobutyrate transaminase 4 Short name= OsGABA-T EC= 2.6.1.-
Expression Region	1-483
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.