



Recombinant Human Differentially expressed in FDCP 8 homolog (DEF8)

Product Code	CSB-BP765060HU
Abbreviation	DEF8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q6ZN54
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MAILSLRAPG PWQAMQVWAD RTLLTPHTGV TSQVLGVAAA VMTPLPGGHA AGRTREARWD AMEYDEKLAR FRQAHLNPFN KQSGPRQHEQ GPGEEVPDVT PEEALPELPP GEPEFRCPER VMDLGLSEDH FSRPVGLFLA SDVQQLRQAI EECKQVILEL PEQSEKQKDA VVRLIHLRLK LQELKDPNED EPNIRVLLEH RFYKEKSKSV KQTCDKCINTI IWGLIQTWYT CTGCYYRCHS KCLNLISKPC VSSKVSHQAE YELNICPETG LDSQDYRCAE CRAPISLRGV PSEARQCDYT GQYYCSHCHW ND LAVIPARV VHNWDFEPRK VSRCSMRYLA LMVSRPVLRL REINPLLFSY VEELVEIRKL RQDILLMKPY FITCREAMEA RLLLQLQDRQ HFVENDEMYS VQDLLDVHAG RLGCSLTEIH TLFAKHIKLD CERCQAKGFV CELCREGDVL FPFDSHTSVC ADCSAVFHRD CYYDNSTTCP KCARLSLRKQ SLFQEPGPDV EA
Source	Baculovirus
Target Names	DEF8
Protein Names	Recommended name: Differentially expressed in FDCP 8 homolog Short name= DEF-8
Expression Region	1-512
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.