



Recombinant *Geobacter sulfurreducens* 50S ribosomal protein L18 (rpIR)

Product Code	CSB-EP765504GBK
Abbreviation	rpIR
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q749A3
Product Type	Recombinant Protein
Immunogen Species	<i>Geobacter sulfurreducens</i> (strain ATCC 51573 / DSM 12127 / PCA)
Purity	>85% (SDS-PAGE)
Sequence	MSPLAQKKVA HLKRKTRVRK KIRGTTERPR LNVFKSAQHI YAQIIDDTTG TTLVSASTVQ KEASESLNHC GNVEAAKRVG ALIAKRALEK NIATVVFDRN GFLYHGRIKA LADSARENGL TF
Source	E.coli
Target Names	rpIR
Protein Names	Recommended name: 50S ribosomal protein L18
Expression Region	1-122
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.