



Recombinant *Ashbya gossypii* Ribosomal N-lysine methyltransferase 5 (RKM5)

Product Code	CSB-YP765591DOT
Abbreviation	RKM5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q750W3
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MLQLCPLDED TLYTHVYERY VELDRHAETL AQDLGIHASD AETLVVDIAP AQPTKSSDTY SLTVSQSLAS LRSSTVNNNS TTGYVLWSGT PFFLCWLLYA PSAAPLRDGG RVPVTDSAAQ FLQLPPLFSA PARPVCVVEL GSGAAGVAAI VLANYVDRYL VSDQKAILKP LRANLLANIS EVSRRTVCSK QTPELSSNRR TPARCELELI ALDWERIATV PAALRPTDAA HVHVLALDVV YNDFLIPPLL TAIKRLLRWY ADEHAVKATA YVLVHLRAQD ILQTFLEHAI IDLRLRCYYM DEERLRSSRF ALYYVTL
Source	Yeast
Target Names	RKM5
Protein Names	Recommended name: Ribosomal N-lysine methyltransferase 5 EC= 2.1.1.-
Expression Region	1-317
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.