



Recombinant *Ashbya gossypii* Protein ATC1/LIC4 (ATC1)

Product Code	CSB-BP765592DOT
Abbreviation	ATC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q751A7
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSVEGKNVDS YSGLTASVQH EHTHPVHNRG VALETTTQHA FITVASQLKP DLKGEDEEHN LELELHKFLS QHNHGPPSSG TTGTTLAGVP DFSASHNISH RDEMAEAIEK AMAEIACPEL WNNDHSSAGA TVGGDVHLDL LTMDGILGNE EPFVQVHEGA LGSVDMDGQR YPHGRNSLSS IVGSKAIIQQ SMDIHPTAQQ GARSRQSTEE VGMHICATPN GKRKNAGVAS SQATVKRSRL FYETISPEL SPLSDNSDFL KAMDTKGSV QNVPISTSVT AGALGHGAKN IKGYQNASKI NMA SNIPKEK LQVHSTPPVV SILPEKLTQE FTMQQVMETK RRIINTHKLI LNFNFLKESY TRSCSELKRT VFCLKESSECH RARLAKENEQ LKRLVIELNE RMKNPSK
Source	Baculovirus
Target Names	ATC1
Protein Names	Recommended name: Protein ATC1/LIC4
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.