



Recombinant *Ashbya gossypii* Serine/threonine-protein phosphatase 2A activator 1 (RRD1)

Product Code	CSB-YP765604DOT
Abbreviation	RRD1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q753L9
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MEFSVEGTQF SPPVKRIFDT AGTHDFQKSL TMYRLQSHLE RYLKLVQGQK IPKSSQNRAV VRFVCILERL DALMDETPPR TGSARRFGDL ACRDWHDRMQ GELDGLLETL LPEAARRSAA ELRYYLGSFAF GSRERLDYGT GHELAF LAVV VALDMLGLWT EDKFTGEDML YVWARYYALV HRLILTYNLE PAGSHGVWGL DDHLHLAYIL GASQWAQDRN VPMQPSDILD PKAVARYSET NLYCNSIAFL LRVKTGHFAQ HSPMLHDIAQ TVPTWSKVTT GLIKMYRVEV LNKFPVVQHF WFGTGFFPWV DMAHGMSLPN YEAPSETSEK PAAGTAHTTT TTMPPPRMTA NCGYGPLGRL VTPRR
Source	Yeast
Target Names	RRD1
Protein Names	Recommended name: Serine/threonine-protein phosphatase 2A activator 1 EC=5.2.1.8 Alternative name(s): Peptidyl-prolyl cis-trans isomerase PTPA-1 Short name= PPIase PTPA-1 Short name= Rotamase PTPA-1 Phosphotyrosyl phospho
Expression Region	1-365
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.