



Recombinant *Ashbya gossypii* Cytoplasmic tRNA 2-thiolation protein 1 (NCS6)

Product Code	CSB-EP765618DOT-B
Abbreviation	NCS6
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q755T1
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MYIPPSDPRR PAKVTAGQLC ELCHARKALV KRPKNLQKVC KLCFFHFVET EIHNTIMENK LFRGERVAV GASGGKDSTV LAYILKLLNE RHDYGLEIVL LSIDEGIVGY RDDSLATVKR NQEYGLPLK IVSYKDLYDW TMDEIVACAG MRNSCTYCGV FRRQALDRGA AMLDIHHVVT GHNADDMAET VLMNILRGDV ARLEKSTAIL TQSSGSPVKR SKPFKYAYQK EIVLYAHYKK LDYFSTECYS APEAFRG TAR ELMKNLEAVR PSCIIDIHS GEALRLRPRP KKRAPPPGSV EIRADGSASL FRNEGFVDGN RCERCGYLSS NRICKACMLL EGGLEKNRARV QIADDTSTEG AARLTRTLEK LHF
Source	E.coli
Target Names	NCS6
Protein Names	Recommended name: Cytoplasmic tRNA 2-thiolation protein 1 EC= 2.7.7.- Alternative name(s): Cytoplasmic tRNA adenylyltransferase 1
Expression Region	1-373
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.