



Recombinant *Ashbya gossypii* Adenylyltransferase and sulfurtransferase UBA4 (UBA4)

Product Code	CSB-MP765623DOT
Abbreviation	UBA4
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q756K6
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MTGESLDGSL HALTIELDAL RRENANLKQQ LKEKDGACGE LPMSLEEFQR YGRQMIVGET GGLSGQVKLR SARVLIVGAG GLGCPALQYL AGAGIGHLGI VDNDVVEESN LHRQPLHDT S KVGLLKDSA KEALSRLNPY CSIKTYPVRL SYANAFEIFP SWDLILDCTD SPMSRYLISD VAVNLGKTVV SGSGLGTEGQ LSIYNFENKG PCYRCFYPIP PRPGSVVSCQ SGGVLGPCIG VLGIMMAVEA LKILFGIYTL ENFKPFLMQY SGFPYQTLRM FKMRNRKQGC LCCGDNPTIT KSTIESGHIK YEAFCGAINY DVLSKDERLS ASEFEANYWS QKERGFVCLD VRPRLHYEIS HLPGTYNMTV KELDEMEGSI EELQKHIPVI TPDLDIVVLC RYGNSRLAT RILKDKFKLR NVRDVKGGYF AYIDEINPSL PKY
Source	Mammalian cell
Target Names	UBA4
Protein Names	Recommended name: Adenylyltransferase and sulfurtransferase UBA4 Alternative name(s): Ubiquitin-like protein activator 4 Including the following 2 domains: Adenylyltransferase UBA4 EC= 2.7.7.- Sulfurtransferase UBA4 EC= 2
Expression Region	1-443
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.