



Recombinant *Ashbya gossypii* Phosphoglycerate kinase (PGK1)

Product Code	CSB-YP765629DOT
Abbreviation	PGK1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q757Q0
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSLSSKLTVK DLSLAGKRVF IRVDFNVPLD GKTITSNQRI VAALPTIKYV LEQGPKAVVL ASHLGRPNGE RNEKYS LAPV AA ELEKLLGQ KVNFLDDCVG EHVTA AVNGA AAGSVF LLEN LRFHIEEEGS RKVDG EKVKA SAEDVQKFRQ GLMSLADV YV NDAFGTAHRA HSSMVGFELP ERAGGFLLSR ELEYFSKALE NPTRPFLAIL GGAKVADKIQ LIDNLLDKVD SIVIGGGMAF TFKKVLENME IGNSIYDKAG AEIVPKLAEK AKKNGVKIVL PVDFVIGDDF SQDANTKIVS ASEGIPSGWE GLDCGPESRK LFSETIASAK TIVWNGPPGV FEIPKFSEGT QAMLAAAVKA SEAGSTVIIG GGDTATVAKK YGVVEKISHV STGGGASLEL LEGKDLPGVT FLSSKQ
Source	Yeast
Target Names	PGK1
Protein Names	Recommended name: Phosphoglycerate kinase EC= 2.7.2.3
Expression Region	1-416
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.