



Recombinant *Ashbya gossypii* Pyruvate kinase (PYK1), partial

Product Code	CSB-MP765637DOT
Abbreviation	PYK1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q759A9
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MESRLGWLTR LNYDTGSALR RTSIIGTIGP KTNNPETLVE LRKAGLNIVR MNFSHGSY EY HQSVIDNARK SEELYPGRPL AIALDTKGPE IRTGTTVDNV DYPIPPNHEM IFTTDDQYAK SCDDKVMFVD YKNITKVIEK GRVIYVDDGV LSFEVLEVVD DKTLKVKSLN AGKICSHKGV NLPGTDVDLP ALSEKDKSDL RFGVKNGVHM VFASFIRTAG DILTIREILG EDGKDIIA KIQNQQGVNN FDEILKVTDG VMVARGDLGI EIPAPQVFAV QKMLIAKCNL AGKPVICATQ MLESMTYNPR PTRAEVSDVG NAVLDGADCV MLSGETAKGI YPINA VKTMA ETALIAEQAI PYVPTYDDLRLNTPKPTPTS ETIAAAAVAA TFEQKAKAIV VLSTSGDTTR LVSKYKPNVP IVMVTRNPR ARYCHLHRGV FPFVYEQEAA DEWTTDDVEAR LNFVGEKSKE LGILSEGDTI VTIQGFASGV GHSNTRLRVLT V
Source	Mammalian cell
Target Names	PYK1
Protein Names	Recommended name: Pyruvate kinase Short name= PK EC= 2.7.1.40
Expression Region	1-501
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.