



Recombinant *Ashbya gossypii* Nuclear distribution protein PAC1 (PAC1)

Product Code	CSB-EP765642DOT-B
Abbreviation	PAC1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q759U7
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MSQHHILPQH QRADLNRSIC DYVQRCGAE SLVVGLRQLF GLSAEECDAT RNGRDPDLLL KKWNSIIRLH RKILDLEQKC QQLTEELEAV PTEAYGKAGR GVSHVLWTPR SNPTFQLDVG ASVTDIKLHP SLPIAFLATD QGKVVAYDLF NRSMPHSTQ AHMKGVTSLA VMEAENGSLI ISTTSKDLQC KIWSFTDNAA FQLLRTLSSH EHIVSQSCFL RSGSDLLLAT CSRDLTVKVV DTKSGWCICKS FQPHNQWVRT LELHGDYVIT GSNDATIRLS HWPSGNGLSM AVMHEFPIER VLIIPMRANT AQKTEADDQ VELEYRKYDP DYSPLGFKYC ISCSRDNLIV LWKIPLPKFI PHRPPQPNNL QTNFEKVHVF KGHTSWVRDI KVRGRHLFSC SDDRSIRCWD LVTGQCLKVW PEASHGFINC LSVSDVSND KLLRELFLSG SIDGKCNVFM R
Source	E.coli
Target Names	PAC1
Protein Names	Recommended name: Nuclear distribution protein PAC1 Alternative name(s): Lissencephaly-1 homolog Short name= LIS-1 nudF homolog
Expression Region	1-461
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.