



Recombinant *Ashbya gossypii* tRNA (guanine (9)-N1)-methyltransferase (TRM10)

Product Code	CSB-EP765644DOT
Abbreviation	TRM10
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q75A17
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	MTPETNNDDET LSRPKPRAAL PPVPEGMSKS QWKKQWKKEQ FELNKPLYAK IRKEKKQKAR EQRRLQKA LEENGGEIPE ELRRTPRVNV NQKDSGIKVI IDCAFDELMN EKEIVSLSTQ ITRAYSANKR ENHFADVKVT SFNKRLKERF DCGLKGANYD AWKHFEFTDE SALPTTNAVY LTADTDETL TLEPGTTYIV GGIVDKNRHK ALCYNKAKEL GIPTRRLPIG EYIKLCGRKV LTTTHVIQIM LRYFDNHDWK EAFESVLPAR KLAELADHAQ ESNSSSPAEE QDAQDI
Source	<i>E.coli</i>
Target Names	TRM10
Protein Names	Recommended name: tRNA (guanine(9)-N1)-methyltransferase EC= 2.1.1.221 Alternative name(s): tRNA methyltransferase 10 tRNA(m1G9)-methyltransferase Short name= tRNA(m1G9)MTase
Expression Region	1-296
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.