



Recombinant *Ashbya gossypii* Elongin-A (ELA1)

Product Code	CSB-MP765678DOT
Abbreviation	ELA1
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q75E67
Product Type	Recombinant Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Purity	>85% (SDS-PAGE)
Sequence	METGHPPSLQ ELCVVVLMRH HAQLEDIGHM PYALVRRVLL KMSAEQLLRL ERASPALLLE DEEAWQQLLK KDFPTNVHEV WVARKKDVYE YYLESVRELG PQVLEDRALV RSLLRRAVCK DPVLFKYRVP SRQLYQRYVH EDARRQELSA ARLRQSTLQL QQEKQKTRIT QLEDPLYLER ELNAGRQRAA PQRSSIYLKS FKELRQRQAH FRSGGYDPTQ RRIVRVQTPS AQPAASLHPP PQSAPMAGSP QKPAPPEQPL AASPTPEAPG PLAQRRRPSA RRQPPPASPA KRRSALFSSP ALSTLLPQQD DTSDSDTGGRP PKKPRVYIHT PR
Source	Mammalian cell
Target Names	ELA1
Protein Names	Recommended name: Elongin-A
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>