



Recombinant Human Centrosomal protein of 68 kDa (CEP68)

| | |
|--------------------------|---|
| Product Code | CSB-MP765741HU |
| Abbreviation | CEP68 |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q76N32 |
| Product Type | Recombinant Protein |
| Immunogen Species | Homo sapiens (Human) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MALGEEKAEA EASEDTKAQS YGRGSCRERE LDIPGPMSE QPPRLEAEGG LISPVWGAEG IPAPTCWIGT DPGGPSRAHQ PQASDANREP VAERSEPALS GLPPATMGSG DLLLSGESQV ETKLSSSEE FPQTLSPRT TTICSGHDAD TEDDPSLADL PQALDLSQQP HSSGLSCLSQ WKSVLSPGSA AQPSSCSISA SSTGSSLQGH QERAEPGGG LAKVSSSLEP VVPQEPSSVV GLGPRPQWSP QPVFSGGDAS GLGRRRLSFQ AEYWACVLPD SLPPSPDRHS PLWNPKEYE DLLDYTYPLR PGPQLPKHLD SRVPADPVLQ DSGVDLDSFS VSPASTLKSP TNVSPNCPA EATALPFGSP REPSLKQWPS RVPQKQGGMG LASWSQLAST PRAPGSRDAR WERREPALRG AKDRLTIGKH LDMGSPQLRT RDRGWPSRP EREKRTSQA RRPTCTESRW KSEEEVESDD EYLALPARLT QVSSLVSYLG SISTLVTLPT GDIKGQSPLE VSDSDGPASF PSSSSQSQLP PGAALQGSGD PEGQNPCFLR SFVRAHDSAG EGS LGSSQAL GVSSGLLKTR PSLPARLDLW PFSDPDVEGQ LPRKGGEQK ESLVQCCKTF CCQLEELICW LYNVADVTDH GTAARSNLTS LKSSLQLYRQ FKKDIDEHQS LTESVLQKGE ILLQCLENT PVLEDVLGRI AKQSGELESH ADRLYDSILA SLDMLAGCTL IPDKKMAAM EHPCEGV |
| Source | Mammalian cell |
| Target Names | CEP68 |
| Protein Names | Recommended name: Centrosomal protein of 68 kDa Short name= Cep68 |
| Expression Region | 1-757 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final |



concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.