



Recombinant Uncharacterized 8.8 kDa protein in frd-Gp32 intergenic region (frd.3)

Product Code	CSB-MP765748BEAA
Abbreviation	frd.3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q76VG8
Product Type	Recombinant Protein
Immunogen Species	Enterobacteria phage LZ4 (Bacteriophage LZ4)
Purity	>85% (SDS-PAGE)
Sequence	MAKVDIDIVD FEYIEEIIRN RYPELSITSI HDDPNYCNXX IVIEGPLDL TRFMANEYCD GMDSEDAEFY MGLIEQ
Source	Mammalian cell
Target Names	frd.3
Protein Names	Recommended name: Uncharacterized 8.8 kDa protein in frd-Gp32 intergenic region
Expression Region	1-76
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.