



Recombinant Chicken Transcription factor Maf (MAF)

Product Code	CSB-MP765766CH
Abbreviation	MAF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q789F3
Product Type	Recombinant Protein
Immunogen Species	Gallus gallus (Chicken)
Purity	≥85% (SDS-PAGE)
Sequence	MASELAMSGS DLPTSPLAME YVNDFDLMKF EVKKEPVETD RIISQCGRLI AGGSLSTPM STPCSSVPPS PSFSAPSPGS GTDQKTHLED YYWMTGYPQQ LNPEALGFSP EDAVEALINS SHHPLPGAFA GYARGQQQLAA AAGGSVPAEE MGSAAAVVSA VIAAAAAQGG APHYHHHHHH PHHGGGGGGG GHPHGAAPGS APPSSASSA AGSGGGGGGG GGGAGGLHHP HHGGGGGGGG LHFDDRFSDE QLVTMSVREL NRQLRGVSKE EVIRLKQKRR TLKNRGYAQS CRFKRVQQRH VLESEKNQLL QQVEHLKQEI SRLVRERDAY KEKYEKLVSN GFRENGSSSD NPSSPEFFM
Source	Mammalian cell
Target Names	MAF
Protein Names	Recommended name: Transcription factor Maf Alternative name(s): V-maf musculoaponeurotic fibrosarcoma oncogene homolog
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.