



Recombinant Mycobacterium tuberculosis

Methionine aminopeptidase 1 (map-1)

Product Code	CSB-MP765906MVZ
Abbreviation	map-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P9WK20
Product Type	Recombinant Protein
Immunogen Species	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
Purity	≥85% (SDS-PAGE)
Sequence	MRPLARLRGR RVVPQRSAGE LDAMAAAGAV VAAALRAIRA AAPGTSSLS LDEIAESVIR ESGATPSFLG YHGYPASICA SINDRVVHGI PSTAEVLAPG DLVSIDCGAV LDGWHGDAAI TFGVGALSDA DEALSEATRE SLQAGIAAMV VGNRLTDVAH AIETGTRAAE LRYGRSFGIV AGYGGHGIGR QMHMDPFLPN EGAPGRGPLL AAGSVLAIEP MLTLGTTKTV VLDDKWTVT ADGSRAAHWE HTVAVTDDGP RILTLG
Source	Mammalian cell
Target Names	map-1
Protein Names	Recommended name: Methionine aminopeptidase 1 Short name= MAP EC= 3.4.11.18 Alternative name(s): Peptidase M
Expression Region	1-266
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.