



Recombinant Danio rerio Pre-mRNA-splicing regulator WTAP (wtap)

Product Code	CSB-MP766455DIL
Abbreviation	wtap
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7SXL7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	MTNEEPLPKK VRLSESDMKT LTREELCTRW KQHETYVQML EAKYADLNSN DVTGLKESEE KKKQQQESA RRENILVMRL ATKEQEMQEC TTQIQYLKQV QQPSAAQLRS SMVDPAINLF FLKMKAEELEQ TKDKLEQAQN ELSAWKFTPD SQTGKKLMAK CRMLIQENQE LGRQLSQGRI AQLEAELALQ KKYSEELKSS QDELNDFIIQ LDEEVEGMQS TILVLQQQLR ETRQQLSQMN QTQGTSSGAG PSRTSPSTAS EPSTQSEPAN ASSSNVGKDC GRVSNGPSNG NSSQRGASGS SLYREASSAD EDYPPSPSVS SPTHDGISKL SNHSEDAVSQ RGGEGYVTQL SAGYESVDSP TGSETSVTQH SNTDTSNADS HEAAAVPKGS RTAGTRHSTQ NGLDSSAAAV ATNTSNASAG SVL
Source	Mammalian cell
Target Names	wtap
Protein Names	Recommended name: Pre-mRNA-splicing regulator WTAP Alternative name(s): Female-lethal(2)D homolog WT1-associated protein Wilms tumor 1-associating protein
Expression Region	1-423
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.