



Recombinant *Oryza sativa* subsp. *japonica* Cysteine protease ATG4B (ATG4B)

Product Code	CSB-MP767123OFG
Abbreviation	ATG4B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XPW8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MTSLPDRGVS SSSSDPLCEG NIAPCSSSSE QKEDCSLKQS KTSILSCVFN SPFNIFEAHQ DSSANKSPKS SSGSYDWSRV LRRIVCSGSM WRFLGTSKVL TSSDVWFLGK CYKLSSEESS SDDSSESGHA TFLEDFSSRI WITYRRGFDA ISDSKYTSDV NWGCMVRSSQ MLVAQALIFH HLGRSWRRPL EKPYNPEYIG ILHMFQDSEA CAFSIHLLQ AGNSYGLAAG SWVGPYAMCR AWQTLVRTNR EQHEVVDGNE SFPMALYVVS GDEDGERGGA PVVICIDVAAQ LCCDFNKGQS TWSPILLVLP LVLGLDKINP RYIPLLKETF TFPQSLGILG GKPGTSTYIA GVQDDRALYL DPHEVQMAVD IAADNIEADT SSYHCSTVRD LALDLIDPSL AIGFYCRDKD DFDDFCSRAT ELVDKANGAP LFTVVQSVQP SKQMYNQDDV LGISGDGNIN VEDLDASGET GEEEWQIL
Source	Mammalian cell
Target Names	ATG4B
Protein Names	Recommended name: Cysteine protease ATG4B EC= 3.4.22.- Alternative name(s): Autophagy-related protein 4 homolog B
Expression Region	1-478
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.