



Recombinant *Oryza sativa* subsp. japonica Serine/threonine-protein kinase SAPK7 (SAPK7)

Product Code	CSB-BP767124OFG
Abbreviation	SAPK7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XQP4
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MERYELLKDI GAGNFGVARL MRNKETKELV AMKYIPRGLK IDENVAREII NHRSLRHPNI IRFKEVVVTP THLAIVMEYA AGGELFDRIC NAGRFSEDEA RYFFQQLICG VSYCHFMIQC HRDLKLENTL LDGSPAPRLK ICDFGYSKSS LLHSPKSTV GTPAYIAPEV LSRREYDGKT ADVWSCGVTL YVMLVGAYPF EDPDDPKNFR KTIGRIMSIQ YKIPEYVHVS QDCRQLLSRI FVANPAKRIT IREIRNHPWF LKNLPRELTE AAQAMYYKGD NSAPTYSVQS VEEIMKIVVE ARTPPRSSTP VAGFGWQEED EQEDNSKKPE EEQEEEEEDAE DEYDKQVKQV HASGEFQLS
Source	Baculovirus
Target Names	SAPK7
Protein Names	Recommended name: Serine/threonine-protein kinase SAPK7 EC= 2.7.11.1 Alternative name(s): Osmotic stress/abscisic acid-activated protein kinase 7
Expression Region	1-359
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.