



Recombinant *Oryza sativa* subsp. *japonica* Putative zinc finger CCCH domain-containing protein 29 (Os04g0487500, LOC_Os04g41060)

Product Code	CSB-EP767127OFG
Abbreviation	Os04g0487500, LOC_Os04g41060
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XSB2
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MTESTKIKYC RADSGSGSKP LSAETKTRLM ADGKKIDEEK KDGTAVRRVL FWFAVLVIRF ALLMAMHYH DVTSSRCDGA RAGGARTGRS ALSRVRGRRR GGGTPDTRPY CDSGTA CLDF RKGCCNHGDA CEFAHGVFDN SSADLA
Source	E.coli
Target Names	Os04g0487550
Protein Names	Recommended name: Putative zinc finger CCCH domain-containing protein 29 Short name= OsC3H29
Expression Region	1-146
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.