



Recombinant *Oryza sativa* subsp. *japonica* Expansin-B5 (EXPB5)

Product Code	CSB-YP767129OFG
Abbreviation	EXPB5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XT39
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	YEQNYT AGRRSTMSLG RGYGWSSGGA TWYGGPQGDG SEGGACGYQS AVGQRPFSSM IAAGGPSLFK NGKGCSCYQ IKCTGNRACS GRPVTVVITD SCPGGVCLNE AAHFDMGTA FGAMANRGMG DRLRSAGVLK IQYKRVPCRF AMNVAFKVDA GSNPYYLAIL VQYANGDGD LAAVHIMEARG GGGWKAMQQS WGATWRLNSN TGKPLSPFFS IRLTSGSGKV LVANNVIPSG WQAGLTYRST VNYAA
Source	Yeast
Target Names	EXPB5
Protein Names	Recommended name: Expansin-B5 Alternative name(s): Beta-expansin-5 OsEXPB5 OsaEXPb1.19
Expression Region	25-275
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.