



Recombinant *Oryza sativa* subsp. japonica Homeobox-leucine zipper protein HOX22 (HOX22)

Product Code	CSB-MP767133OFG
Abbreviation	HOX22
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7XUJ5
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MNGRTQLASW ARIAMDRGDH HHLQQQHQL MPPPAPVVPP QLCMPAMMAD EQYMDLGGGG AAAAPGRGGA GERKRRFTEE QIRSLESMFH AHHAKLEPRE KAELARELGL QPRQVAIWFQ NKRARWRSKQ LEHDYAALRS KYDALHSRVE SLKQEKLALT VQLHELRLERL REREERSGNG GAATTAASSS SCNGSGSEEV DDDDDKRNAA AGCLDLEPPE SCVLGGATCA TPADVSVESD QCDDQLDYDE GLFPESFCAT PELWEPWPLV EWNAVA
Source	Mammalian cell
Target Names	HOX22
Protein Names	Recommended name: Homeobox-leucine zipper protein HOX22 Alternative name(s): HD-ZIP protein HOX22 Homeodomain transcription factor HOX22 OsHox22
Expression Region	1-276
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.