



# Recombinant Halocynthia roretzi Histone-lysine N-methyltransferase SETD7 (setd7)

<b>Product Code</b>	CSB-EP767185HTQ-B
<b>Abbreviation</b>	setd7
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7Z0G7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Halocynthia roretzi (Sea squirt) (Cynthia roretzi)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDSSDDEIAC DEG DYKGAKD DNDLPHGLGK VKFSSGDEFI GAFEHGIKCG PGKFHFFDDS TLEGN YVDGE LHGIGIYTND DGSITKSTYC EGVMEGPSWE YDPEGNITFR GQYSEGVR CG LCFYFPDGG SLIGNVNASG DLSADNIAYI YDRTTALIG SFEEGDMITA KEANVTITGE KGEEISFPTV NSISPDPVYR LDVSTPHVIS TRPLVPDPYE SELVYAAPSK IPNAGEGLYA KCDVDQDTVM AFYNGVRLKQ DEVENRDWSQ NSNTISLTDD IAIDVPEEYV STDNYCASLG HKVNHSFDPN CRYDIYQHPR FGFIKCVRTI RGVSEGDEL T VHYTYEHNDG NKTREAEAPE WYKSQLKVFG VDRPAEILEN MDEDYC
<b>Source</b>	E.coli
<b>Target Names</b>	setd7
<b>Protein Names</b>	Recommended name: Histone-lysine N-methyltransferase SETD7 EC= 2.1.1.43 Alternative name(s): SET domain-containing protein 7
<b>Expression Region</b>	1-386
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.