



Recombinant Human ATP synthase subunit g 2, mitochondrial (ATP5L2)

Product Code	CSB-EP767202HU-B
Abbreviation	ATP5L2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z4Y8
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MAPFVRNLVE KTPALVNAAV TYLKPRLAAF WYYTTVELVP PTPAEIPRAI QSLKKIVSSA QTGSFKQLTV KEALLNGLVA TEVSTWFYVR EITGKRGIIG
Source	E.coli
Target Names	ATP5L2
Protein Names	Recommended name: ATP synthase subunit g 2, mitochondrial Short name= ATPase subunit g 2
Expression Region	1-100
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.