



# Recombinant Human F-box only protein 33 (FBXO33)

<b>Product Code</b>	CSB-EP767209HU-B
<b>Abbreviation</b>	FBXO33
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q7Z6M2
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MLLFLSVPQP RPPGARTRAG AARVARWRRL RLQQLRRLRG LLRVLRGRPG AGSRRRGRMA LCGQAAGAAS LPSELIVHIF SFLPAPDRLR ASASC SHWRE CLFYPALWPQ LRICLRVSPA EQPRLEFLMR KCGWVRELR VEFAAENYLS GGGPGDGGGA DTGTGGEEVE ALQLSARWLE VLRTYLELVL CVLVSIRNNR NLQKFSLFGD ISVLQQQGS LNTYLSKVDP DGKKIKQIQQ LFEEILSNR QLKWLSCGFM LEIVTPTSLS SLSNAVANTM EHL SLLDNNI PGNSTLITAV ELERFVNLHS LALDFCDFTA EMARVLTDSN HVPLQRLSLL VHNVSVMHKS LDNMPNDEHW KALSRKSTSF RYIIMAFDIK SEDMLKILKP SIPLERIHFD SYITCVSGAI VDLISRQYDK FLTHFILMND VIDTSGFPDL SDNRNEDPLV LLAWRCKLS LLAIHGYTVW AHNLI AIARL RGS DLKVLEV TEESIDFDQG ELADQDVDPV HNLIEQVSLG LGQPWHAVMD IESLSVFTEP NRHFYREMQS FSEDI
<b>Source</b>	E.coli
<b>Target Names</b>	FBXO33
<b>Protein Names</b>	Recommended name: F-box only protein 33
<b>Expression Region</b>	1-555
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.