



Recombinant Human F-box only protein 33 (FBXO33)

Product Code	CSB-MP767209HU
Abbreviation	FBXO33
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z6M2
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	MLLFLSVPQP RPPGARTRAG AARVARWRRL RLQQLRRLRG LLRVLRGRPG AGSRRRGRMA LCGQAAGAAS LPSELIVHIF SFLPAPDRLR ASASC SHWRE CLFY PALWPQ LRICLRVSPA EQPRLEFLMR KCGWVRELR VEFAAENYLS GGPGDGGGA DTGTGGEEVE ALQLSARWLE VLRTYLELVL CVLVSIRNNR NLQKFSLFGD ISVLQQQGS LNTYLSKVDP DGKKIKQIQQ LFEEILSNR QLKWLS CGFM LEIVTPTSLS SLSNAVANTM EHL SLLDNNI PGNSTLITAV ELERFVNLHS LALDFCDFTA EMARVLTDSN HVPLQRLSLL VHNVSVMHKS LDNMPNDEHW KALSRKSTSF RVYIMAFDIK SEDMLKILKP SIPLERIHFD SYITCVSGAI VDLISRQYDK FLTHFILMND VIDTSGFPDL SDNRNEDPLV LLAWRCKLS LLAIHGYTVW AHNLI AIARL RGS DLKVLEV TEESIDFDQG ELADQDVDPV HNLIEQVSLG LGQPWHAVMD IESLSVFTEP NRHFYREMQS FSEDI
Source	Mammalian cell
Target Names	FBXO33
Protein Names	Recommended name: F-box only protein 33
Expression Region	1-555
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.