



Recombinant Schizosaccharomyces pombe Uncharacterized ABC transporter ATP-binding protein C20G4.01 (SPAC20G4.01)

Product Code	CSB-EP767218SXV
Abbreviation	SPAC20G4.01
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7Z991
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	>85% (SDS-PAGE)
Sequence	MEVTVSNLSY TFSPKQPLSL DHVTLDLPKG SRTLLVGANG AGKSTLLKLL SGKSLAKAGH ISVGGKDPFR ESSSAFVYLG TEWVNNPVIH RDMSVARLIA SVGGDKFAER RDFLISILDI DLRWRMHAVS DGERRRVQLC MGLLRPFVFL LLDEVTVDLV VLARADLLNF LQEETEVRNA TIVYATHIFD GLAEWPTHLV HLSLGRIVDY GPISKFGALM TRSSTGNSAL LETCLEWLKE DKKNRGTREE EKRSTWDDVQ ESLKVGTDVF TDYFKISRAK
Source	E.coli
Target Names	SPAC20G4.01
Protein Names	Recommended name: Uncharacterized ABC transporter ATP-binding protein C20G4.01
Expression Region	1-280
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.