



Recombinant Danio rerio Acylglycerol kinase, mitochondrial (agk)

Product Code	CSB-MP767241DIL
Abbreviation	agk
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q7ZW00
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	>85% (SDS-PAGE)
Sequence	NHCDSV LRREACIEAR AFGQQLIGPQ EILKKATVIL NPAACKGKAN QLFEKNAAPI LHLAGVEVKI VKTDYEGQAK KLMELMEQTD MLIAGGDGT LQEVITGLLR RADEEIFSKT PIGFIPLGSS NSLSQSLHLV SDNKVQHITS ATLSILKGET VPLDVLQIKG EKEQPVFALL GLRWGAFRDV ATSISKYWYL GPLKTRAAHW FSSLKQWPQV HVASLSYLAP VPRPPDLPDE IPQRPNLLYR IYRRLQNYWN PPLEEPPKEP EPEQWESKDV STELTLVLTH NKNPVKRRVD DSMLISLDPE SLTAGQFITE GTNKSVDPME PIENAVQIEA SAARLELPEE GAGFYDIDNQ EYEAMSVEVR LLPRKLRFFC SAERREQLAE AQ
Source	Mammalian cell
Target Names	agk
Protein Names	Recommended name: Acylglycerol kinase, mitochondrial EC= 2.7.1.107 EC= 2.7.1.94 Alternative name(s): Multiple substrate lipid kinase Short name= MuLK Short name= Multi-substrate lipid kinase
Expression Region	35-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.