



Recombinant Danio rerio Protein Wnt-10b (wnt10b)

Product Code	CSB-EP767270DIL-B
Abbreviation	wnt10b
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q801F7
Product Type	Recombinant Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Purity	≥85% (SDS-PAGE)
Sequence	ND ILGLKVAGEP VLTPNAVCLR LAGLTKKQMR LCVRSPDVTA SALQGIQVAI HECQHQLRDQ RWNCSSLENH GKLPHQSAIL NRGFRESAFS LSLLAAGVVH SVASACSLGK LRGCGCEAKR RLDDDKIRLK LTQLQLQTFQ RSGVSLAGAG ENTPELSSLH GSLPANLHSS HPMSLLKPLP DEVTMLQDTW EWGGCSDIR FGVRFSDWL DSRGSPRDIH ARTRIHNNRV GRQVVTDNMR RKCKCHGTSG SCQFKTCWYV SPEFRLVGSL LREKFLTAIF INSQNKNGV FNSRTGGSTG SDPLRGQRRR SISRELVYFE KSPDFCDREP AVDSLGTQGR ICNKSSPGMD GCGSLCCGRG HNILKQARSE RCHCRFWCC YVLCEECKVT EWWNVCK
Source	E.coli
Target Names	wnt10b
Protein Names	Recommended name: Protein Wnt-10b
Expression Region	29-427
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.