



Recombinant Oreochromis mossambicus Natriuretic peptides A (nppa)

Product Code	CSB-MP767282OES
Abbreviation	nppa
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q805E9
Product Type	Recombinant Protein
Immunogen Species	Oreochromis mossambicus (Mozambique tilapia) (Tilapia mossambica)
Purity	>85% (SDS-PAGE)
Sequence	KRTSGCF GARMDRIGNA SGLGCNSGRG
Source	Mammalian cell
Target Names	nppa
Protein Names	Recommended name: Natriuretic peptides A Alternative name(s): Prepronatriodilatin Cleaved into the following 2 chains: 1. Atrial natriuretic factor Short name= 2. ANF Alternative name(s): Atrial natriuretic peptide Short name=
Expression Region	114-140
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Cytoplasmic domain
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.