



Recombinant Mouse Vacuolar protein sorting-associated protein 37D (Vps37d)

Product Code	CSB-MP767372MO
Abbreviation	Vps37d
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q810I0
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	MYRARAARAG PEPGSPGRFG ILSTGQLRDL LQDEPKLDRI VRLSRKFQGL QLERDACLAS NYALAKENLA LRPRLEMGR T ALAIKYQELR EVAENCADKL QRLEKSMHRW SPQCALGWLQ AELEEAQEQA EVQMEQLLLG EQSLEAF LPA FQRGRAL AHL RRTQAEKLQE VLRRRERSAQ PAPTTAAAAA AAATAMDPPK PFPAAVLPT GAARGPPPAV PRSLPPLDSR PVPPVKGSPG CPF GPAPLLS PRPSQPEPPH R
Source	Mammalian cell
Target Names	Vps37d
Protein Names	Recommended name: Vacuolar protein sorting-associated protein 37D Alternative name(s): ESCRT-I complex subunit VPS37D Williams-Beuren syndrome region protein 24 homolog
Expression Region	1-261
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.