



Recombinant Chlamydomonas reinhardtii UPF0301 protein CCA_00630 (CCA_00630)

Product Code	CSB-YP767590DSL
Abbreviation	CCA_00630
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q822P9
Product Type	Recombinant Protein
Immunogen Species	Chlamydomonas reinhardtii (strain GPIC)
Purity	>85% (SDS-PAGE)
Sequence	MAKIPYAILE KGSLLLASPD TDQGVFARVS ILLCEHSLNG SFGLILNKTL GLEISDDIFT FDKVSNNNIR FCMGGPLQAN QMMLLHSCSE ISEQTLEICP SVYLGGLSF LQEIAASESG PTINLCFGYS GWQAGQLEKE FLEGNWFLAP ASYEYVFSDN PDNLWSRILK DLGGKYASLS TVPENLLL
Source	Yeast
Target Names	CCA_00630
Protein Names	Recommended name: UPF0301 protein CCA_00630
Expression Region	1-189
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.