



Recombinant Chlamydomonas reinhardtii Triosephosphate isomerase (tpiA)

Product Code	CSB-YP767599DSL
Abbreviation	tpiA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q823U7
Product Type	Recombinant Protein
Immunogen Species	Chlamydomonas reinhardtii (strain GPIC)
Purity	≥85% (SDS-PAGE)
Sequence	MERKSYVFGN WKMHKTAKEA KDFLSIFCPF VKEISPASLV GIAPAF TTLS ACCESIKEMD SSIWLGAQNV HQDSSGAFTG EVSLPMLQEF HVNFVLLGHS ECRHIFHEED ATIALKVGAA ARSGVVPVLC IGETLEAREN GTTKEVLSNQ LMLGLAQLPE TAPVIIAYEP VWAIGTGKVA STADVQEVHA FCREVLSRIF SKEKSEIISI LYGGSVKADN AEGFARCPDV DGLLVGGASL DPQGFANVLR NFNL
Source	Yeast
Target Names	tpiA
Protein Names	Recommended name: Triosephosphate isomerase Short name= TIM EC= 5.3.1.1 Alternative name(s): Triose-phosphate isomerase
Expression Region	1-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.