



Recombinant Chlamydomonas reinhardtii UPF0109 protein CCA_00022 (CCA_00022)

Product Code	CSB-EP767611DSL
Abbreviation	CCA_00022
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q824W6
Product Type	Recombinant Protein
Immunogen Species	Chlamydomonas reinhardtii (strain GPIC)
Purity	>85% (SDS-PAGE)
Sequence	MKDFLAYIIK NLVDRPEEVH IKEVQGTHTI IYELTVAKPD IGKIIGKEGR TIKAIRTLV SVASRNNVKV SLEIMEDK
Source	E.coli
Target Names	CCA_00022
Protein Names	Recommended name: UPF0109 protein CCA_00022
Expression Region	1-78
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.