



Recombinant UPF0678 fatty acid-binding protein-like protein SAV_4030 (SAV_4030)

Product Code	CSB-YP767656SNN
Abbreviation	SAV_4030
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q82G68
Product Type	Recombinant Protein
Immunogen Species	Streptomyces avermitilis (strain ATCC 31267 / DSM 46492 / JCM 5070 / NBRC 14893 / NCIMB 12804 / NRRL 8165 / MA-4680)
Purity	>85% (SDS-PAGE)
Sequence	MIEIPSDLHK DLVPLAFLLG NWAGAGVHDF PGSEKCNFGQ EVSFTHDGRD FLEYHSHTWV LDADGNKVRP LETESGFWRV DAARKVEAVM TRDDGVVEIW YGELAAKKPQ IDLVTDVAR TAAAGPYSGG KRLYGYVKSD LMWVGEKQTP EVELRPYMSA HLKKVTPED VERWAKALPD DMPDDGIAFF K
Source	Yeast
Target Names	SAV_4030
Protein Names	Recommended name: UPF0678 fatty acid-binding protein-like protein SAV_4030
Expression Region	1-191
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.