



Recombinant Nitrosomonas europaea UPF0042 nucleotide-binding protein NE1849 (NE1849)

Product Code	CSB-YP767692NHH
Abbreviation	NE1849
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q82TN5
Product Type	Recombinant Protein
Immunogen Species	Nitrosomonas europaea (strain ATCC 19718 / CIP 103999 / KCTC 2705 / NBRC 14298)
Purity	>85% (SDS-PAGE)
Sequence	MQVIIIISGLS GSGKSIALKV LEDSGYYCVD NLPASLLVVL INHLQTQQHA YVAVAIDMRS GENITVLPWQ LKMIDKSIQI KFIFLEARTE TLMQRFSETR RRHPLSDKNI TLEEAIRRER EALATLTGLG HHIDTSSLRP NVLRAFIKDF IADSRSPSQL TLLFQSFGYK HGIPLDADLV FDIRCLPNPF YDPQLKELTG HDPEVIRFME SQPDASKMLR DISSFLGTWL PAYIRDNRAY LTVAIGCTGG QHRSVYFAEK LALHFHDSAHLVLRHRGLAE YKPHYARR
Source	Yeast
Target Names	NE1849
Protein Names	Recommended name: UPF0042 nucleotide-binding protein NE1849
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.