



Recombinant Nitrosomonas europaea Putative ribosome biogenesis GTPase RsgA (rsgA)

Product Code	CSB-YP767730NHH
Abbreviation	rsgA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q82Y12
Product Type	Recombinant Protein
Immunogen Species	Nitrosomonas europaea (strain ATCC 19718 / CIP 103999 / KCTC 2705 / NBRC 14298)
Purity	>85% (SDS-PAGE)
Sequence	MSVGADRKKT GTSLTGKVVA AYGRHFEVEV AGGTIYSCVV RGKKKGVVCG DEVEILPATG DQGIIETTLR RTSLFYRSEI FREKLIANA TQLVFLAVV PSCNLELLDR CLVAAESQGI RPLILLNKID LIGQDEQRQA VAHHLMFYRE LGYPVLEISA KISVQPLIPL LSGQTSLLAG QSGVVGKSTLL NALVPRAQQA TAEISDALDS GRHTTTHVRL FHFADASSII DSPGFQEFGL QQLDEASLAR GFIEFRPFLG QCKFRDCRHI AEPGCKLLLA AQEGMLNSRR IACYHKLKVG LKKSHPWMET NKRK
Source	Yeast
Target Names	rsgA
Protein Names	Recommended name: Putative ribosome biogenesis GTPase RsgA EC= 3.6.1.-
Expression Region	1-314
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.