



# Recombinant Pantothenate kinase (coaA)

<b>Product Code</b>	CSB-MP767936SZB
<b>Abbreviation</b>	coaA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q83PC5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Shigella flexneri
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSIKEQTLMT PYLQFDRNQW AALRDSVPMT LSEDEIARLK GINEDLSLEE VAEIYLPLSR LLNFYISSNL RRQAVLEQFL GTNGQRIPYI ISIAGSVAVG KSTTARVLQA LLSRWPEHRR VELITTDGFL HPNQVLKERG LMKKKGFPEP YDMHRLVKFV SDLKSGVPNV TAPVYSHLIY DVIPDGDKT VQPDILILEG LNVLQSGMDY PHDPHHVFVS DFVDFSIYVD APEDLLQTWY INRFLKFREG AFTDPDSYFH NYAKLTKEEA IKTAMTLWKE INWLNKQNI LPTRERASLI LTKSANHVVE EVRLRK
<b>Source</b>	Mammalian cell
<b>Target Names</b>	coaA
<b>Protein Names</b>	Recommended name: Pantothenate kinase EC= 2.7.1.33 Alternative name(s): Pantothenic acid kinase
<b>Expression Region</b>	1-316
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.