



Recombinant 1- (5-phosphoribosyl)-5-[(5-phosphoribosylamino)methylideneamino]imidazole-4-carboxamide isomerase (hisA)

| | |
|--------------------------|---|
| Product Code | CSB-MP767954SZB |
| Abbreviation | hisA |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q83R04 |
| Product Type | Recombinant Protein |
| Immunogen Species | Shigella flexneri |
| Purity | >85% (SDS-PAGE) |
| Sequence | MIIPALDLID GTVVRLHQGD YGKQRDYGND PLPRLQDYAA QGAEVLHLVD LTGAKDPAKR QIPLIKTLVA GVNVPVQVGG GVRSEEDVAA LLEAGVARVV VGSTAVKSPD VVKGWFERFG ADALVLALDV RIDEQGNKQV AVSGWQENSG VSLEQLVETY LPVGLKHVLC TDISRDGTLA GSNVSLYEEV CARYPQVAFQ SSGGIGDIND VAALRGTGVR GVIVGRALLE GKFTVKEAIA CWQNA |
| Source | Mammalian cell |
| Target Names | hisA |
| Protein Names | Recommended name: 1-(5-phosphoribosyl)-5-[(5-phosphoribosylamino)methylideneamino]imidazole-4-carboxamide isomerase EC= 5.3.1.16 Alternative name(s): Phosphoribosylformimino-5-aminoimidazole carboxamide ribotide isomerase |
| Expression Region | 1-245 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life |



of lyophilized form is 12 months at -20°C/-80°C.