



# Recombinant Ferrochelatase (hemH)

<b>Product Code</b>	CSB-MP767967SZB
<b>Abbreviation</b>	hemH
<b>Storage</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
<b>Uniprot No.</b>	Q83SE5
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Shigella flexneri
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRQTKTGILL ANLGPDPAPT PEAVKRYLKQ FLSDRRVVDV S RLLWWPLL R GVILPLRSPR VAKLYASVWM EDGSPLMVYS RQQQQALAGR LPDTPVALGM SYGSPSLESS VDELLAEHVD HIVVLPYPQ FSCSTVGAVW DELARILARK RSIPGISFIR DYADNHDIYIN ALANSVRASF AKHGEPDLLL LSYHGIPQRY ADEGDDYPQR CRTTTRELAS ALGMVPEKVM MTFQSRFGRE PWLMPYTDDET LKMLGEKGVG HIQVMCPGFA ADCLETLEEI AEQNREVFLG AGGKKYEYIP ALNATPEHIE MMANLVAAAYR
<b>Source</b>	Mammalian cell
<b>Target Names</b>	hemH
<b>Protein Names</b>	Recommended name: Ferrochelatase EC= 4.99.1.1 Alternative name(s): Heme synthase Protoheme ferro-lyase
<b>Expression Region</b>	1-320
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>