



Recombinant Lipid-A-disaccharide synthase (lpxB)

Product Code	CSB-YP767968SZB
Abbreviation	lpxB
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q83SL3
Product Type	Recombinant Protein
Immunogen Species	Shigella flexneri
Purity	>85% (SDS-PAGE)
Sequence	<pre>MTDQRPLTIA LVAGETSGDI LGAGLIRALK ERVPNARFVG VAGPRMQAEG CEAWYEMEEL AVMGIVEVLG RLRLLHIRA DLTKRFGELK PDVFGIDAP DFNITLEGNL KKQGIKTIHY VSPSVWAWRQ KRVFKIGRAT DLVLAFLPFE KAFYDKYNVP CRFIGHTMAD AMPLDPDKNG ARDVLGIPYD AHCLALLPGS RGAEVEMLSA DFLKTAQLLR QTYPDLEIVV PLVNAKRREQ FERIKAEVAP DLSVHLLDGM GREAMVASDA ALLASGTAAL ECMLAKCPMV VGYRMKPFTE WLAKRLVKTD YVSLPNLLAG RELVKELLQE ECEPQKLAAL LLPLLANGKT SHAMHDTFRE LHQQIRCND EQAAQAVLEL AQ</pre>
Source	Yeast
Target Names	lpxB
Protein Names	Recommended name: Lipid-A-disaccharide synthase EC= 2.4.1.182
Expression Region	1-382
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>