



Recombinant Regulatory protein E2 (E2)

Product Code	CSB-EP767977HOM
Abbreviation	E2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q84294
Product Type	Recombinant Protein
Immunogen Species	Human papillomavirus type 6a
Purity	≥85% (SDS-PAGE)
Sequence	MEAIKRLDA CQEQLLELYE ENSTDNLKHV LHWKCMRHES VLLYKAKQMG LSHIGMQVVP PLKVSEAKGH NAIEMQMHLE SLLKTEYSME PWTLQETSYE MWQTPPKRCF KKRKGTVEVK FDGCANNTMD YVWVTDVYVQ DTDSWVKVHS MVDAGIYYT CGQFKTYYN FVKEAEKYGS TKQWEVCYGS TVICSPASVS STTQEVSIP STTYTPAQTS TPVSSSTQED AVQTPPRKRA RGVQQSPCNA LCVAHIGPVD SGNHNLITNN HDQHQRNNNS NSSATPIVQF QGESNCLKCF RYRLNDKHRH LFDLISSTWH WASPKAPHKH AIVTVTYHSE EQRQQFLNVV KIPPTIRHKL GFMSLHLL
Source	E.coli
Target Names	E2
Protein Names	Recommended name: Regulatory protein E2
Expression Region	1-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.