



# Recombinant Beta-glucanase (licB)

<b>Product Code</b>	CSB-BP767985CMO
<b>Abbreviation</b>	licB
<b>Storage</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
<b>Uniprot No.</b>	Q84C00
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Clostridium thermocellum (Ruminiclostridium thermocellum)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	EAA TVVNTPFVAV FSNFDSSQWE KADWANGSVF NCVWKPSQVT FSNGKMILT L DREYGGSSYPY KSGEYRSKSF FGYGYEVRM KAAKNVGI SFFTYTGPSD NNPWDEIDIE FLGKDTTKAQ FNWYKNEVGG NEYLHNLGFD ASQDFHTYGF EWRPDYIDFY VDGKKVYRGT RNIPVTPGKI MMNLWPGKGV DEWLGRYDGR TPLQAEYEYV KYYPNGVPQD NPTPTPTIAP STPTNPNLPL KGDVNGDGHV NSSDYSLFKR YLLRVIDRFP VGDQSVADV N RDGRIDSTDL TMLKRYLIRA IPSL
<b>Source</b>	Baculovirus
<b>Target Names</b>	licB
<b>Protein Names</b>	Recommended name: Beta-glucanase EC= 3.2.1.73 Alternative name(s): 1,3-1,4-beta-D-glucan 4-glucanohydrolase Endo-beta-1,3-1,4 glucanase Laminarinase Lichenase
<b>Expression Region</b>	28-334
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>