



Recombinant Beta-glucanase (licB)

Product Code	CSB-EP767985CMO
Abbreviation	licB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q84C00
Product Type	Recombinant Protein
Immunogen Species	Clostridium thermocellum (Ruminiclostridium thermocellum)
Purity	≥85% (SDS-PAGE)
Sequence	EAA TVVNTPFVAV FSNFDSSQWE KADWANGSVF NCVWKPSQVT FSNGKMILT DREYGGSSYPY KSGEYRSKSF FGYGYEVRM KAAKNVGI SFFTYTGPSD NNPWDEIDIE FLGKDTTKAQ FNWYKNEVGG NEYLHNLGFD ASQDFHTYGF EWRPDYIDFY VDGKKVYRGT RNIPVTPGKI MMNLWPGKGV DEWLGRYDGR TPLQAEYEYV KYYPNGVPQD NPTPTPTIAP STPTNP NLPL KGDVNGDGHV NSSDYSLFKR YLLRVIDRFP VGDQSVADV N RDGRIDSTDL TMLKRYLIRA IPSL
Source	E.coli
Target Names	licB
Protein Names	Recommended name: Beta-glucanase EC= 3.2.1.73 Alternative name(s): 1,3-1,4-beta-D-glucan 4-glucanohydro lase Endo-beta-1,3-1,4 glucanase Laminarinase Lichenase
Expression Region	28-334
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.