



Recombinant *Oryza sativa* subsp. *japonica* 14-3-3-like protein GF14-A (GF14A)

Product Code	CSB-MP767992OFG
Abbreviation	GF14A
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q84J55
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAAAAGGGTR EEMVYMAKLA EQAERYEEMV EFMEKVVTTAA AAGGGGELTV EERNLLSVAY KNVIGARRAS WRIVSSIEQK EEGRGAAGHA AAARSYRARV EAELSNICAG ILRLDLERLV PAAAVDAKV FYLKMKGDYH RYLAEFKTGA ERKDAADATL AAYQAAQDIA MKELSPHPI RLGLALNFSV FYYEILNSPD RACTLAKQAF DEASELDTL GEESYKDSTL IMQLLRDNL LWTSDMQDDG GDEMRRDATKP EDEH
Source	Mammalian cell
Target Names	GF14A
Protein Names	Recommended name: 14-3-3-like protein GF14-A Alternative name(s): G-box factor 14-3-3 homolog A
Expression Region	1-264
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.