



Recombinant Arabidopsis thaliana Cytokinin riboside 5'-monophosphate phosphoribohydrolase LOG8 (LOG8)

Product Code	CSB-MP768003DOA
Abbreviation	LOG8
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q84MC2
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MEDNQRSRFR KICVFCGSHS GHREVFSDAA IELGNELVKR KIDLVIYGGGS VGLMGLISRR VYEGGLHVLG IIPKALMPIE ISGETVGDVR VVADMHERKA AMAQEAEAFI ALPGGYGTME ELLEMITWSQ LGIHKKTVGL LNVDGYNNL LALFDTGVVEE GFIKPGARNI VVSAPTAKEL MEKMEEYTPS HMHVASHESW KVEELGDYPG QENKPQ
Source	Mammalian cell
Target Names	LOG8
Protein Names	Recommended name: Cytokinin riboside 5'-monophosphate phosphoribohydrolase LOG8 EC= 3.2.2.n1 Alternative name(s): LOG family protein At5g11950 Protein LONELY GUY 8
Expression Region	1-216
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.