



Recombinant Arabidopsis thaliana F-box protein At3g12350 (At3g12350)

Product Code	CSB-EP768028DOA
Abbreviation	At3g12350
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q84WW1
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MATFALPFCE IPEDLQLRIL SLLTPAEISS FACTSKRFAS LCQEDGKIWH VMCDQRWGKK TKIQKWANGQ IPYRHLYKTL KGLENLIGFW RLCGRANPAA SSPPLVFFDW GSSFILGSRV LSTGDDTYEV KKTPFLLMGI SPEGRSENFL DLVGGNLRVS DDDLKELEAS NNLVSVDVNF MGNGHIMVEE NRCFTNNRR EDQRSSGDES DDLISSPNFS EMYTQLANKT SPGGDRRRQK RKEKERQASR TKWEPEHFLK VADCSPTPTR PLQGLWKGFC EGSIELYLVK YDEVGGIICR KVEDLSLSRY TPPVFWTPKH TFIRSPFSAE EELLNSRIH ISPFAEVHEN VVSGMLYMK SCDLVLPGEP GNGIGFLRGE GRVWLYDNGT FGFGFLRDQF IIDLKRVALE DGCLADEVEA CM
Source	E.coli
Target Names	At3g12350
Protein Names	Recommended name: F-box protein At3g12350
Expression Region	1-422
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.