



Recombinant *Oryza sativa* subsp. *japonica* Putative B3 domain-containing protein LOC_Os07g12820 (LOC_Os07g12820)

Product Code	CSB-YP768031OFG
Abbreviation	LOC_Os07g12820
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q84YL3
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MHSPTFSMVK IKTASQDYLP IPVAVTKASR LKHGRTLKLM TAHGLKIRVK VAEARDKLYM TIGWKEFIQE AGLKMGESKS VVFRTLSKSR LNVIIFNKEG YSRCPIPDKA AKALINNQSS SAPSFSTKST APRHPSFTNV EANTKRIVKD MCCYNKRMKL SSEVKNYVRD IAQFLDYSSK FYIVTMKNIH EVRQGDGCYN LKAGYGQISH KIIFKEIYSA TSTASQQDTS GKNLHFSLPQ DESVAGNSTS PVNAMFMPCT SNLLLL
Source	Yeast
Target Names	Os07g0232100
Protein Names	Recommended name: Putative B3 domain-containing protein LOC_Os07g12820
Expression Region	1-267
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.